

# PEAK ELLIPSE HEAVY DUTY

ROLLATOR & WHEELCHAIR - 9346376000765

\$770



KG

205kg WEIGHT  
CAPACITY



FOLDS FOR  
EASY TRANSPORT



FULLY ADJUSTABLE  
SUSPENSION



# PEAK ELLIPSE HEAVY DUTY

ROLLATOR & WHEELCHAIR - 9346376000765

The Peak Ellipse Heavy Duty rollator features adjustable handles, a large, supportive backrest and a deep seat, and folds for easy transport. The Ellipse Heavy Duty has a 205kg weight capacity, solid puncture-proof tyres and the ability to be used as a rollator or wheelchair. The Peak Ellipse Heavy Duty rollator is perfect for travelling and built for comfort.

## SPECIFICATIONS

<b>Adjustable Backrest</b>	Yes	<b>Number of components</b>	1
<b>Adjustable Seat Height</b>	No	<b>Seat Dimensions</b>	584 x 220
<b>Adjustable Handle Height</b>	915 - 1040mm	<b>Seat height from floor</b>	600mm
<b>Dimensions (L x W x H)</b>	495 x 610 x 1040mm	<b>Seat padding</b>	Firm
<b>When folded (L x W x H)</b>	495 x 226 x 915mm	<b>Total height</b>	915 - 1040 (adjustable)
<b>Distance between armrests</b>	584mm	<b>Tyres</b>	Solid Puncture Proof
<b>Drive Wheel</b>	204mm	<b>Weight Capacity</b>	205kg
<b>Front Wheel</b>	254mm	<b>Total Weight</b>	9.3kg
<b>Heaviest part</b>	11kg	<b>Main Material</b>	Aluminium